


Monitoring Gastric Emptying via a Simple Non-invasive Breath Test

^{13}C breath tests offer a non-invasive, non-radioactive alternative to many other diagnostic techniques. Gastroparesis (delayed gastric emptying) is often diagnosed via scintigraphy - a costly, cumbersome procedure that must be performed in specialized outpatient centers and exposes the patient to radiation-emitting compounds. Alternatively, the rate of emptying can be much more easily determined by administering the patient with a ^{13}C labeled substrate, and monitoring the expelled $^{13}\text{CO}_2$ in the patient's breath on a Sercon ABCA.

The Cairn ^{13}C -Spirulina Gastric Emptying Breath Test (GEBT) is the only FDA-approved breath test used in the diagnosis of gastroparesis. The kit contains a specially formulated meal of dehydrated scrambled eggs containing pharmaceutical-grade spirulina, enriched in ^{13}C . The Sercon ABCA is a high performance, high sample throughput isotope ratio mass spectrometer, the accepted gold standard for breath testing since 1999 and offers the highest accuracy and performance available for ^{13}C breath testing.



Breath samples are collected from the patient, before and after eating Cairn's GEBT meal, and analysed via isotope ratio mass spectrometry on a Sercon ABCA. After administration of the meal, the labeled by-products of the ^{13}C -Spirulina (proteins, carbohydrates and fats) are absorbed and metabolized, giving rise to $^{13}\text{CO}_2$ which is expelled in the breath. The rate of $^{13}\text{CO}_2$ excretion is proportional to the rate of gastric emptying. This simple non-invasive test is an effective means of measuring gastric emptying, and has been validated against scintigraphy.

Cairn Diagnostics GEBT using Sercon's ABCA is FDA approved and both the GEBT kit and Sercon's ABCA are CE marked medical devices. Cairn Diagnostics analyses breath samples in a CLIA certified lab.

^{13}C	Sample	Specification (Standard Deviation)
Ref Gas Injection Precision	100% CO_2 (n=5)	$\leq 0.05\%$
Ref gas in tubes precision	10ml 5% CO_2 (n=5)	$\leq 0.05\%$
Breath Precision	Breath in Exetainer (n=5)	$\leq 0.1\%$
Linearity	10 to 20ml 5% CO_2 in Exetainer	$\leq 0.3\%$ change



All enquiries regarding testing services at Cairn should be addressed to info@cairndx.com



info@sercongroup.com
www.sercongroup.com